INSTRUKTION EASYLIFT SHOULDER

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Organisation number
556086-9413

INSTRUCTION **EASYLIFT SHOULDER**

EasyLift an equipment which is designed to ease the burden especially on the users fingers, hands and arms when lifting/carry. The advantage is rapidly accruing easing of weight from your fingers, hands and arms to the shoulders. Handling of the stretcher is similar to lifting and handling without EasyLift equipment, be careful to follow the instructions for safe lifting and unloading of heavy duty and carry over longer distances.

1. Make sure to dress EasyLift equipment with harness and buckles on the front.

2. Make sure to dress with the cross on the back located in the middle of upper back.


3. Adjust length to fit your body by placing the middle finger in the finger pocket at the bottom of carry straps on each side. Adjust the straps so that you have slightly bented arms.

4. Grasp the handle a bit longer towards your hand wrist than normal.

5. When you clench your hand may your fingers slip out from the finger pocket, this does not affect the use, finger pocket is only to guide the strap to the start position.

6. On the friction strap, at the bottom there is a Velcro sewn, fold up the strap and attach it. You can now wear EasyLift over or under your jacket out of sight. If you get too much strap out, would it be an advantage to fold it into the carabines to be positioned not getting in the way.

Distribution:  **FERNO**
Norden

Address
Box 2127
SE-291 02 Kristianstad
SWEDEN
info@germa.se
www.germa.se

Street
Industrigatan 54-56

Phone
+46(0)44123030

Telefax
+46(0)44103179

Organisation number
556086-9413