DIRECTIONS FOR USE

VACUUM SPLINTS

English: Directions for use. **Ferno Vacuum Splints – Ferno AS - Orange Line**
Deutsch: Bedienungsanleitung. **Ferno Vakuumschienen – Ferno AS - Orange Line**
Français: Mode d'emploi. **Attelles à dépression Ferno AS - Série orange**
Español: Instrucciones de uso. **Férulas de Vacío de Ferno AS - Línea naranja**
Svenska: Bruksanvisning. **Ferno AS Vakuumsplintar - Orange modell**
Italiano: Istruzioni per l'uso. **Stecche a Depressione Ferno AS - Linea arancio**
Português: Instruções de utilização. **Talas de Vácuo Ferno AS - Linha laranja**
Dansk: Brugsanvisning. **Ferno AS Vakuum Splints - Orange type**
Nederlands: Gebruiksaanwijzing. **Ferno AS Vacuüm Spalken zijn - Orange Line**
Suomi: Käyttöohje. **Ferno AS Vacuum Splint - Oranssi linja**
Norsk: Brukerveiledning. **Ferno AS Vacuum Splint - Oransje serie**

Distributor:

Ferno Washington Inc.
70 Weil Way
Wilmington, OH 45177
USA

Phone: +1 937 382 1451
Fax: +1 937 382 6569

E-mail: info@ferno.com

www.ferno.com

Manufacturer:

AB Germa
Industrigatan 54-56
SE-291 36 Kristianstad
Sweden

Phone: +46 (0)44 12 30 30
Fax: +46 (0)44 10 31 79

E-mail: info@germa.se

www.germa.se

English:	Directions for use	Page	3
Deutsch:	Bedienungsanleitung	Seite	7
Français:	Mode d'emploi	Page	11
Español:	Instrucciones de uso	Página	15
Svenska:	Bruksanvisning	Sida	19
Italiano:	Istruzioni per l'uso	Pagina	23
Português:	Instruções de utilização	Página	27
Dansk:	Brugsanvisning	Side	31
Nederlands:	Gebruiksaanwijzing	Bladzijde	35
Suomi:	Käyttöohje	Sivu	39
Norsk:	Brukerveiledning	Side	43

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1. Warning and caution statements

CAUTION

- Check regularly that the vacuum pump works correctly and that the air can be evacuated from the Vacuum Splint so as to make it rigid.
- Avoid handling the Vacuum Splint close to sharp objects.
- Avoid handling the Vacuum Splint close to an open fire.
- Check regularly that the vacuum pump, the Vacuum Splint and the straps are not damaged or display signs of wear and tear. Repair any damages and change defective parts immediately.
- Check regularly that the Vacuum Splint has no cracks.

2. Intended use

The **Ferno Vacuum Splints** are intended to stabilise injured patients during transport. They are suited for fixation of patients with hand, arm, leg, shoulder and hip injuries. The Vacuum Splints are available as full leg-, lower leg-, arm-, forearm- and wrist splints.

3. Specifications

The Ferno Vacuum Splints are in conformity with ISO 6940:1984
The Ferno Vacuum Splints are in conformity with Council Directive 93/42/EEC concerning Medical Device

Temperature range: -30 °C to +70 °C

Parts/Materials

Cover material: PVC coated polyester fabric (flame retardant)
Inner net: Polyester
Granules: Polystyrene granules (flame retardant)

Dimensions	Length (cm)	Width (cm)	Height (cm)	Weight (kg)	Capacity for hand pump
Ferno Full leg AS 150	127	98/50	3	1,30	-450hPa at c. 40 strokes
Ferno Leg AS 140	112	36/70	3	1,00	-450hPa at c. 25 strokes
Ferno Arm AS 120	75	48/63	3	0,65	-450hPa at c. 14 strokes
Ferno Forearm AS 100	55	36	3	0,40	-450hPa at c. 9 strokes

4. Description

The Vacuum Splint consists of an outer cover made of hard-wearing PVC coated fabric with an inner net of polyester filled with polystyrene granule. The inner net is welded inside the splint to secure granules stay in the splint.

When the air is evacuated, the granule is moulded into a tight substance making the Vacuum Splint rigid and unable to shape. Before and while the air is evacuated the Vacuum Splint may be moulded according to the shape of the patient's limb in order to obtain optimum support of the patient.

The Vacuum Splint is equipped with straps with Velcro fastenings and for Full Leg splint - 2 extension straps.

5. Operating instructions

Full Leg Splint – see 1, Leg Splint – see 2, Arm Splint – see 3, Forearm Splint – see 4, Wrist Splint – see 5.

1. Put the splint on an even surface near the patient.
2. Mould the splint so that it is even or, if necessary, in an appropriate state for the emergency.
3. Close the valve by turning it to the right (clockwise).
4. Place the leg/lower leg/arm/forearm/wrist in the medically desired position on the splint.
5. Mould the splint around the leg/arm so that the fracture becomes stable and the unit is comfortable for the patient.
6. Attach the Velcro fasteners and tighten them.
7. Lift the foot part of the splint to a supporting angle of the foot and fasten the strap crosswise to the leg.
8. Connect the pump. Evacuate the air, and continue to mould the splint. Tighten the Velcro fastenings once more.
9. The patient is ready to be transported.

For patients with short legs the foot part should be folded.

By using the extension straps it is possible to support fractures in hip or shoulder.

In special cases, the Full Leg Splint can be used as a vacuum mattress for a small child. For this you need the optional head strap.

In special cases, the Leg Splint can be used as a vacuum mattress for a baby. For this, you need the optional head strap.

6. Cleaning and disinfecting

After each use the Vacuum Splint should be wiped with a moist cloth and dried. If dirty the splint may be washed with soapy water and rinsed with clean water.

Allow the Vacuum Splint to dry completely before storing it.

If the Vacuum Splint needs disinfection first clean it carefully and let it dry. Wipe the complete splint with a solution containing 70% alcohol or use the disinfectant called LYSETOL AF.

Test the Vacuum Splint at regular intervals to make sure that it has the necessary tightness.

If leaks are found; see section 7 for repair instructions.

7. Maintenance and repair

At a regular basis the Vacuum Splint and the pump must be checked for visible damages.

At the same time the splint must be tested as to tightness. Furthermore the suction force of the vacuum pump must be checked according to Section 8.

Leaks can be found on the Vacuum Splint by visual inspection or by filling the splint with air covering it with soapy water all over the surface and the welds and watching out for air bubbles from the splint.

If leaks have been found on the splint, mark the leak, let any air out, clean the area around the leakage and let it dry.

Find a patch in a suitable size and put glue on the patch as well as on the leakage on the splint. Wait for the glue to become dust dry (approx. 2 minutes) and place the patch on the leakage. Firmly press the whole area of the patch establishing good contact and press out all air under the patch.

The repair will have full strength after 24 hours.

8. Test of the Vacuum Splint and Pump

Straighten out the Vacuum Splint, **close the valve** by turning it clockwise and connect the pump.

Evacuate the air until the Vacuum Splint feels rigid and leave it for about 2 hours; if the splint is as rigid as it was when you started it may be considered tight. Open the valve by turning it counter-clockwise.

If the Vacuum Splint is not airtight, see section 7 for repair or check the valve for leakage. Try to insert a new valve and repeat the test.

If the pump is not emptying the Vacuum Splint check the pump by closing the connection with your fingertip and pump a few times.

A negative pressure should be created in the tubing that can be felt at your fingertip. If the pump does not function do not repair it but return it for service to your dealer or to an Ferno service centre.

9. Storage

Loosen the valve by turning it counter-clockwise and put the Vacuum Splint on an even surface.

When the granules have been spread evenly, fold the Vacuum Splint in its initial position.

Repack the Vacuum Splint so that it fits into your storage area. The Vacuum Splint is often positioned folded in a bag.

10. Ordering information - Ferno Orange Line

005 002 015	Ferno Full Leg Splint AS 150
005 002 014	Ferno Leg Splint AS 140
005 002 013	Ferno Arm Splint AS 120
005 002 012	Ferno Forearm Splint AS 100

11. Spare parts

005 002 017	Repair kit - Orange Line
271 000 716	Multifit Valve
271 000 709	Sealing ring for pump, large size 63,5
271 000 710	Sealing ring for pump, small size 14,6

12. Accessories

271 000 701	Ferno Vacuum foot pump
272 000 701	Ferno Vacuum hand pump
272 000 705	Carrying bag for splints
005 002 118	Ferno AS 190 Splint Set - Orange Line (Forearm AS 100, Arm Splint AS 120, Leg Splint AS 140, hand pump, connector and bag)
272 000 706	Head strap
272 000 705	Extension strap

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