


Lite Model

Article numbers: Retmex Lite Model:	1099602010
Size rescue sheet (L x W):	188x69 cm
Fabric material:	100% Polyamide (complies with DIN 1865-1)
Packaging:	Storage bag attached to Rescue Sheet
Number of required users:	2 - 6 persons
Weight:	+/- 500 gram (complies with DIN 1865-1)
Weight capacity / load limit:	250 kg (complies with DIN 1865-1)
Number of handles:	12 (complies with DIN 1865-1)
Cleaning instructions:	In the washing machine: 40 degrees, dry on low temperature, iron on low temperature, no chemical cleaning, do not bleach. By hand: Use a wipe and disinfectant to clean by hand
Guarantee:	The Rescue Sheets will have a 1-year warranty for material and manufacture faults, which is valid from the date of delivery. The manufacturer (Tetcon) may ask for proof of the purchase date. Therefore, it is recommended to keep the invoice in a safe place. Faults or problems with the use of the Rescue Sheets as a result of an accident, misuse, incorrect use, alterations, vandalism, incorrect placement, incorrect maintenance or administration, as well as overdue normal and routine maintenance, are not covered by the warranty. Tetcon is not liable for (consequential) damage as a result of the use of the Rescue Sheet.
Conformity:	The product are classified in class 1 according to regulation MDR 2017/745 / EU. Applied standards: NEN-EN-ISO 13485: 2016 & NEN-EN-ISO 14971: 2019 EN



Lite Model



Backside

