

# NORSE RESCUE® Non-Rebreathing Mask XL

---



Item No: BR-00113  
 Manufacturer: EMSRUN  
 Certificates: CE, FDA  
 Size: XL  
 Color: Green/Transparent  
 Country of origin: Denmark/China  
 NSN no: N/A  
 Tariff: 90192000  
 Weight: 90 g

NORSE RESCUE® Non-rebreathing Mask supply high oxygen concentration to the patient.

The NORSE RESCUE® Non-Rebreathing Mask comes with a reservoir bag and is used in emergency situations to prevent hypoxemia, also known as low blood oxygen.