

TSL Expander Kit

Expander for Scoop EXL

Codice: 25-0030-002



TSL Expander Kit



TSL Expander is an accessory made from high-strength aluminum alloy that attaches to the Scoop EXL atraumatic stretcher and its standard TSL hook fitting.

Designed for the movement and transport of a large patient in all intra- and extra-hospital environments.

The option of using the Expander Kit, the Scoop EXL atraumatic stretcher, and Easy Lift hangers, also allows you to use the system with lifters.

The kit contains: 2 x Expanders (red), 3 x restraint extensions, and 1 x container bag.



KEY SPECIFICATIONS

- Kit contents can be stowed anywhere.
- High-strength metal restraints
- Assembles in seconds.
- Fully decontaminable.
- Mortuary belts available on request
- Compatible with Easy Lift hangers



Scoop EXL code: 0708002

Code: 0313659-IT-M

Code: PS0765EXL

OPTIONAL ACCESORIES

Description	Code
TSL Expander Kit	25-0030-002
Easy Lift (optional)	PS0765EXL
Mortuary belts (optional)	0313659-IT-M



TECHNICAL SPECIFICATIONS

Height	150 mm
Width	300 mm
Weight (KIT)	1.5 kg
Load capacity	227kg