



VACUUM MATTRESS GERMA EASYFIX PLUS II

Our latest vacuum mattress EasyFix PLUS II is made of flame-retardant PVC-coated fabric. It has new inner bags with chambers for controlling granules, which means that it always has an even distribution of the contents for a comfortable and stable fixation.

EasyFix PLUS II has a unique design that provides increased support for the head, neck and back without the use of a neck collar on the patient. The unique design also facilitates the management of the airways of a patient with limited mobility. Despite this, the packaging volume and weight have decreased by 20% compared to most traditional systems.



Art. Nr	GE23271511001
Certifications	MDR
Color	Rød/sort
Load Capacity	250 kg
Weight	8,7 kg / 11 kg
Height	4,5 cm
Length	208 cm
Width	130/72 cm



KEY FEATURES

- Back - extended area fiberglass shaft - main area
- Reinforced with Poly propylene PP board to reduce wear at transitions to wood board
- Headband - enhanced with anti-slip surface against the patient's forehead
- The position of the Headband - adjusted down towards the back creates a better angle and stabilization of the patient's head