

UNI Transfer Mattress is designed in co-operation with professional paramedic users focusing on ergonomics, patient handling, patient safety and comfort for patient.

UNI Transfer Mattress has been tested by trauma teams to secure efficient handling plus quick, safe and easy transfer of patients including heavy patients.

UNI Transfer Mattress approved according to EN 1865 and EN 1789

Art no: 271 030 003

Features and Benefits:

**Safety**

Easy to adjust safety belts according to standard EN 1865 and EN 1789

**Comfort**

**Extra** comfort by using mix of Sensus memory foam 20 mm and comfort 40 mm V-foam which support patients body weight and deliver a safe and comfort feeling.

**Ergonomics**

Strong handles to support lift, carry or slide of patient on the mattress using 10 strong handles. It is possible to lower the head for upper airway management. Loading and carry capacity 250 kg.

Suitable to be used in combination with Step by Step system.



**Designed to fit on top of most common stretcher designs.**

- **Hygienic** full welded easy to clean with most common disinfectants
- **Sliding material / reducing friction** material under the mattress, easy to slide patient over to e.g. hospital bed.
- \*Fits together with Germa Step by Step sliding kit.



Accessories: Step by Step

Art no: 15 000 000 117



**Specification— Universal Transfer Comfort Mattress**

PU- Polyester fibre flame retardant  
 PP-plate and grass fibre material to reduce friction  
 Foam Sensus — V-shaped (material that deliver extra comfort and side support to the patient)

**Dimensions and weight**

Length: 185 x 50 x 9 cm  
 Weight: 4,7 kg  
 Handles: 10

**Ordering Specification**

UNI Transfer Mattress

**Art. No.**

271030003