



## VACUUM SPLINT, SHORT ARM

Vacuum terminals model Ferno<sup>®</sup> are of high quality in a simple and easy-to-understand design. Vacuum terminal Ferno<sup>®</sup> is one of the world's best-selling Vacuum terminals. The splinters are manufactured in Sweden, but the high turnover makes this splinter very affordable. Ferno<sup>®</sup> vacuum terminals are designed to stabilize and fix fractures in patients during transport

|                |         |
|----------------|---------|
| Art. Nr        | GEAS100 |
| Certifications | MDR     |
| Weight         | 0,75 kg |
| Length         | 55 cm   |
| Width          | 36 cm   |
| Height         | 36 cm   |



### KEY FEATURES

- The model is suitable for fixing patients with arm or wrist injuries.
- The splint forms around the arm and provides optimal and comfortable support.
- Durable pvc fabric makes the splinters durable and easy to clean.
- All seams are fully welded for best hygiene.
- The splinters hold a vacuum for at least 15 hours and are X-ray permeable.
- One-year warranty applies to production defects and quality.