

## Features

- 18 oz. polyester reinforced PVC vinyl and with all seams and zippered attachments being double stitched for durability
- The main zipper is a heavy duty #10 coil zipper with double sliders allowing the top flap to be completely opened for easy loading and unloading
- For use with all models of Traverse Rescue basket stretchers

## Design

The Traverse Rescue Stretcher Cover protects the basket stretcher from the environmental elements, whether you store it on the top of your rescue vehicle or on a job site. It is designed with ample room to allow for the storage of backboards, splints, blankets and other stretcher accessories and comes with 6 reinforced carry handles and two compression straps.

Constructed using 18 oz. polyester reinforced PVC vinyl and with all seams and zippered attachments being double stitched for durability. The main zipper is a heavy duty #10 coil zipper with double sliders allowing the top flap to be completely opened for easy loading and unloading. For use with all models of Traverse Rescue basket stretchers.



Part No.	Description	Width	Length	Depth	Weight
05-0144	Basket Stretcher Cover	66 cm / 26 in	216 cm / 85 in	19 cm / 7.5 in	2.45 kg / 5.5 lb
05-0145	Split Basket Stretcher Cover	66 cm / 26 in	111.76 cm / 44 in	19 cm / 7.5 in	1.45 kg / 3.25 lb