

The Porter Transport System is a fast, simple securing system with 2 never flat wheels and eight position handles, offering an all purpose transport solution. Greater leverage and easier carrying reduces fatigue and chance of injury.

## Features

- Universal design, compatible with any type and style of rigid litter
- Eight position adjustable handles
- Fast ratchet securing system locks basket stretcher with 2 straps and 2 rubberized steel hooks/strap to cinch it down
- Quick release pins on each side of the handle locks them into position with ease; Each pin has a wire lanyard
- Design allows the wheels to be quickly and easily removed for compact storage and transport

## Design

- Compact, easy storage, wheels fit inside frame
- Capable of securing a makeshift litter made from wooden planks or logs
- Two securing straps (2 hooks per strap) are each rated to 4.5 kN (1,000 lbs.) and use a strong and durable Ratchet-Locking device easily operated in all weather conditions
- Hooks on securing straps have a one-inch diameter opening and can be attached to any style litter.



Part No.	Description	Width	Length	Height	Weight	Load Rating
11-2224	The Porter	28 in / 71 cm	37.5 in / 95 cm	19.5 in / 50 cm	41 lbs / 18.5 kg	500 lbs / 227 kg