



# FIRST AID BAG ALS L ROLLER

**Large yet compact respiratory care bag with an interior designed for professional use.**

Thoughtfully designed with ergonomic carry handles and integrated wheels, making it easy to handle even when heavily packed. The integrated trolley with a telescopic handle makes it ideal for travel. The bag has hidden shoulder straps and can also be worn as a backpack.



## KEY FEATURES

- Large exterior pocket on the front of the bag
- Hidden, lockable telescopic handle
- Two ergonomic rubber carry handles
- Hidden pocket with zipper for backpack straps
- Transparent outer pocket
- Reflective strips all around for high visibility
- Heavy-duty zippers with large pulls for easy use with gloves

<b>Art. Nr</b>	VT235
<b>Color</b>	Red
<b>Load Capacity</b>	25 kg
<b>Liquid Capacity</b>	66 L
<b>Material</b>	600D Polyester
<b>Weight</b>	6,2 kg
<b>Depth</b>	31 cm
<b>Height</b>	61 cm
<b>Width</b>	35 cm

